Rosetta-gami 2(DE3)

Rosetta-gami 2(DE3) Chemically Competent Cell 说明书

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D(ara-leu)7697 DlacX74 DphoA Pvull phoR araD139 ahpC galE galK rpsL
(DE3) F' [lac+ lacIq pro] gor522::Tn10 trxB pRARE2 (CamR, StrR, TetR)

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Rosetta-gami 2(DE3) 菌株集合了 Rosetta 2 和 Origami 2 两种菌株的优点:

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Sample Induction Protocol (for reference only)

1. Inoculate a single colony from a freshly streaked plate into 3ml of LB medium containing the appropriate antibiotic for the plasmid and host strain.
2. Incubate with shaking at 200 rpm at 37°C overnight.
3. Inoculate 50 ml of LB medium containing the appropriate antibiotic with 0.5 ml of the overnight culture prepared in step 2 (use the 500 ml triangular flask as the container would be better).
4. Incubate with shaking at 150 rpm at 37°C until the OD 600 reaches 0.5-0.8. (0.6 recommended; about 2.5h).
5. (Optional) Pipet 1ml of the cultures into clean microcentrifuge tubes and place the tubes on ice until needed for gel analysis or storage at -20°C. These will serve as the non-induced control samples.
6. Add IPTG to a final concentration of 1 mM. Optimal time for induction of the target protein may vary from 2-16 hours, depending on the protein.
7. Incubate with shaking at 120 rpm at 37°C for 2-4 hours. To determine the optimal time for induction of the target protein, it is recommended that a time course experiment be performed varying the induction from 2-16 hours.
8. Place the culture on ice for 10 minutes. Harvest cells by centrifugation at 5,000 x g for 10

minutes at 4 °C.

9. Remove the supernatant and store the cell pellet at -20 °C (storage at lower temperatures is also acceptable).

IPTG

Prepare a 1 M solution of IPTG (Isopropyl- β -D-thiogalactoside; Isopropyl- β -D-thiogalactopyranoside) by dissolving 2.38 g of IPTG in dd water and adjust the final volume to 10 ml. Filter sterilize before use.